

Performance of a probe-type graphite calorimeter (Aerrow) in ultra-high dose per pulse electron beams


**Alexandra Bourguin¹, F. Keszti², W. Simon³, A. Schoenfeld³, J. Seuntjens²,
J. Renaud⁴, A. Schüller¹, R.-P. Kapsch¹**

(1) Physikalisch-Technische Bundesanstalt (PTB), Braunschweig, Lower Saxony, DE,

(2) McGill University, Montreal, QC, CA,

(3) Sun Nuclear Corporation, Melbourne, FL,

(4) National Research Council Canada, Ottawa, ON, CA



No, nothing to disclose

- Standard dosimeter for external beam reference dosimetry is the ionization chamber
 - The ultra-high dose rate means that ion recombination is very large > 80% !

- Integrating dosimeters
 - Alanine or radiochromic film have been shown to work
 - BUT not real-time (readout after 24 hours)

- ... may be calorimeter

➤ A probe-type graphite calorimeter



J. Renaud, Med Phys 2018, DOI: 10.1002/mp.12669

- Developed at McGill University
- Provided by Sun Nuclear Corporation
- Designed for clinical measurement
- Quasi-adiabatic mode

➤ Reading system


- Temperature sensor (thermistor) directly read by high-stability DMM Agilent

- Metrological Electron Accelerator Facility (MELAF) at PTB, Germany



- Tests carried out at 20 MeV, 5 Hz PRF, pulse width of 2.5 μ s
- Dose varied between 0.6 Gy and 5.6 Gy per pulse
- Detector response compared to beam current monitor

- Absolute dose to water measurement
 - Absolute dose measurement for a range of dose per pulse
 - Compared to a calibrated Advanced Markus (k_{sat} measured at the PTB)
 - For 4 fixed dose per pulse, tested different number of pulses delivered




➤ Absolute dose to water measurement

- Absolute dose measurement for a range of dose per pulse
- For 4 fixed dose per pulse, tested different number of pulses delivered
- Compared to a calibrated Advanced Markus (k_{sat} measured at the PTB)

➤ Depth dose curve measurement

- Maximum dose rate
- Irradiation time of 2 s (10 pulses delivered)
- Compared to Advanced Markus and diamond detector prototype



Dose to water




$$D_w = c_{gr} \cdot \Delta T \cdot K_{ht} \cdot \left(\frac{D_w}{D_{gr}} \right)_{MC}$$

Dose to water



Specific heat capacity of graphite
Consensus value

$$D_w = c_{gr} \cdot \Delta T \cdot K_{ht} \cdot \left(\frac{D_w}{D_{gr}} \right)_{MC}$$

Dose to water




Specific heat capacity of graphite
 Consensus value

$$D_w = c_{gr} \cdot \Delta T \cdot K_{ht} \cdot \left(\frac{D_w}{D_{gr}} \right)_{MC}$$




Temperature change
Measurement

Dose to water




Specific heat capacity of graphite
Consensus value

$$D_w = c_{gr} \cdot \Delta T \cdot K_{ht} \cdot \left(\frac{D_w}{D_{gr}} \right)_{MC}$$



Temperature change
Measurement

Heat loss correction
Thermal Simulation




Dose to water




Specific heat capacity of graphite
Consensus value

Dose conversion factor
Monte Carlo





$$D_w = c_{gr} \cdot \Delta T \cdot K_{ht} \cdot \left(\frac{D_w}{D_{gr}} \right)_{MC}$$



Temperature change
Measurement

Heat loss correction
Thermal Simulation






➤ Monte Carlo simulation


- EGSnrc
- Dose ratio conversion factor
- Energy mapping for thermal simulation

➤ Thermal simulation


- FEM simulation (COMSOL)
- To calculate k_{ht}



Figures from: J. Renaud, Med Phys 2018, DOI: 10.1002/mp.12669



Results: variation dose per pulse




$$(D_{w/gr})_{MC} = 1.149(1)$$

- Fit
- ◇ Aerrow
- Advanced Markus

UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: variation number of pulses




Dose rate
Gy per pulse

- 0.6
- ◆ 2.2
- ▲ 4.1
- 5.3

UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: variation number of pulses




Dose rate
Gy per pulse

- 0.6
- ◆ 2.2
- ▲ 4.1
- 5.3

UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: variation number of pulses




Dose rate
Gy per pulse

- 0.6
- ◆ 2.2
- ▲ 4.1
- 5.3


UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: depth dose curve




UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: depth dose curve



UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

Results: depth dose curve



UNPUBLISHED DATA-DO NOT COPY OR DISTRIBUTE

➤ Absolute dosimetry

- Calorimeter is showing promising results for NMI and clinical dosimetric use
- Heat lost correction factor about 0.8% (correction factor k_{sat} for ICs about 20-80%)

➤ For relative measurement

- Correction factor not required (in contrast to ICs where a correction function is needed)
- Real time measurement would be possible with specific software



Acknowledgments

Andreas Schlesner
Christoph Makowski
Leon De Prez

<http://uhdpulse-empir.eu/>

This project 18HLT04 UHDpulse has received funding from the EMPIR programme co-financed by the Participating States and from the European Union's Horizon 2020 research and innovation programme.



The EMPIR initiative is co-funded by the European Union's Horizon 2020 research and innovation programme and the EMPIR Participating States



**Physikalisch-Technische Bundesanstalt
Braunschweig**

Bundesallee 100
38116 Braunschweig

Alexandra Bourgouin

Telefon: +49 531 592-6439

E-Mail: alexandra.bourgouin@ptb.de

www.ptb.de

Merci
Thank you
Danke schön