



Schweizerische Eidgenossenschaft  
Confédération suisse  
Confederazione Svizzera  
Confederaziun svizra

Federal Institute of Metrology METAS



## **FRICKE DOSIMETRY AS A PRIMARY STANDARD AND REFERENCE FOR ABSORBED DOSE TO WATER IN ULTRA HIGH PULSE DOSE RATE ELECTRON BEAMS**

F. Frei, M. Trachsel, M. Liechti , P. Peier


# Faculty Disclosure

<input checked="" type="checkbox"/>	No, nothing to disclose
	Yes, please specify:

Company Name	Honoraria/ Expenses	Consulting/ Advisory Board	Funded Research	Royalties/ Patent	Stock Options	Ownership / Equity Position	Employee	Other (please specify)

# Agenda

## 1. Calibration using Fricke solution



**3. Consequences for using Fricke dosimetry in UHD electron beams.**

$$\varepsilon \cdot G$$

measure


or

ICRU90

or

compare

## 2. Determination of radiation chemical yield -> Total absorption technique




to transfer instruments: Alanine, EBT3, ionisation chamber, ...

For 1L Fricke solution @ METAS:

392.14 mg	Fricke salt	$(\text{NH}_4)_2\text{Fe}(\text{II})(\text{SO}_4)_2$
+ 58.44 mg	Sodium chloride	NaCl
+ up to 1L	Aqueous sulfuric acid	$\text{H}_2\text{SO}_4$ , 0.4mol/L

# Calibration using Fricke solution



- Irradiation of the Fricke solution leads to the oxidation of  $\text{Fe}^{2+}$  to  $\text{Fe}^{3+}$ .
- The increase of  $\text{Fe}^{3+}$  is proportional to the absorbed dose.



$$D_W = D_F \cdot f_{F/W} = \dots = \frac{\Delta A \cdot f_{F/W}}{\rho \cdot l \cdot \epsilon(\lambda) \cdot G} = \frac{A_{\text{irradiated}} - A_{\text{unirradiated}}}{\rho \cdot l \cdot \epsilon(\lambda) \cdot G}$$

measured Monte Carlo

Absorbed dose in Fricke solution in Gy

Transfer factor

Change in absorbance  $A_{\text{irradiated}} - A_{\text{unirradiated}}$

Chemical yield in mol/J


Density of Fricke solution in kg/l

Path length through cuvette in cm


Molar extinction coefficient in  $\text{l}/\text{cm}/\text{mol}$

## Fricke dosimeter:


- Closely water-equivalent ( $f_{F/W}$ : 1.005 - 1.008, depending on the electron energy)
- Independent of irradiation energies for energies  $> 2 \text{ MeV}$
- Bag (50  $\mu\text{m}$  thick polyethylene foil) of size 30x40x3 mm<sup>3</sup>



# Determination of radiation chemical yield -> Total absorption



- Special electron beam.
- Glass vessel, closed using a 50µm polyethylene foil (volume: 1L or 2L)
- Temperature stabilised to 25°C using circulating water




$$\varepsilon \cdot G = \frac{\Delta A \cdot f_{F/W}}{\rho \cdot l \cdot D_W} = \frac{\Delta A \cdot f_{F/W}}{[\rho] \cdot [l]} \cdot \frac{m}{[E_{el}] \cdot [N] \cdot f_p}$$

measured  
Monte Carlo


Mass of Fricke solution  
 Electron energy  
 Number of electrons  
 Fraction of incident electron beam energy deposited

- Fraction of incident beam electrons energy deposited in the Fricke solution:  $f_p$ : 92 – 97% (depending on energy and vessel)
- Uncertainty total absorption: aiming for  $\varepsilon \cdot G$  @304nm 1.4% ( $k=2$ )
- Next step: measurement campaign to compare to other primary standards

# Consequences for using Fricke dosimetry in UHD electron beams



- Allows measurements in the «UHD-reference beam»
- At  $\sim 1\text{Gy/pulse}$  (of  $3\mu\text{s}$ ) radical-radical reactions start to reduce the value of radiation yield G.
- Next step: further characterize the Fricke dosimeter at higher instantaneous dose rates to determine correction factors.
- Uncertainty for conventional:  $N_{W,Q} 1.7\%$  ( $k=2$ )



# Acknowledgement

The project 18HLT04 UHDpulse has received funding from the EMPIR programme co-financed by the Participating States and from the European Union's Horizon 2020 research and innovation programme.



The EMPIR initiative is co-funded by the European Union's Horizon 2020 research and innovation programme and the EMPIR Participating States





Schweizerische Eidgenossenschaft  
Confédération suisse  
Confederazione Svizzera  
Confederaziun svizra

Federal Institute of Metrology METAS



Thank you very much for your attention